



# For IV and Epidural Procedures

# 3123

## START CART WITH ARTICULATING ARM with Specialty Package, Four Drawers & Latch Lock

Harloff's Start Cart is designed to reduce the risk of accidental needle sticks to patients and staff while improving the efficiency of the IV or epidural procedure.



Engineered & Manufactured in the U.S.

### Features:

- Heat sealed powder coat steel cabinet & drawers
- Pontoon bumper protects facility walls and cart
- Articulating arm allows clinician to work from sitting or standing position
- Dual lipped working surfaces with clear plastic cover overlay
- Adjustable push handle in 1-inch increments
- Soft grip handles for moving cart from sitting position
- 3 gallon sharps and waste containers make it easy to dispose of supplies
- Ballasted for stability and lower center of gravity
- 16" full extension drawer slides
- Basic latch lock
- Premium 5-inch casters make it easy to roll around the OR and patient rooms
- Removable plastic top

### Benefits:

- Clinicians will no longer have to gather needed equipment from a distant cart, resulting in going back and forth for supplies from the cart to the patient site.
- Clinicians will no longer have to use the patient's bed, stretcher or surgical operating room table as a work surface.
- Clinicians will no longer have to leave exposed needles on top of the patient's bed that may cause injury to the doctor, nurse, medical technician, or housekeeper
- Technicians can prepare supplies to be used in advance by clinicians which will save time and money.



Color shown is Sand. More than 20 additional colors are available.

Dual work surfaces and an articulating arm provide work area to prepare and dispose of IV and Epidural supplies. Cart provides compact storage for point of use supplies.

#### Assembled Cart Dimensions:

- Base cart: 21.50" w x 22.50" d (54.61 cm W x 57.15 cm D)
- Cart height Base Unit: 30" h (76.20 cm H)
- Articulating arm extends height to 45.75" h (116.21 cm H)
- Cart with waste and sharps containers: 31.50" w x 22.50" d (80.01 cm W x 57.15 cm D)
- Adjustable Push Handle: Lowest height 34.75" extends to 39.75" in 1" increments (Metric equivalent: Lowest height 88.27 cm extends to 100.97 cm in 2.54 cm increments)
- Cart net weight: 116 pounds (52.62 kg)
- Approximate shipping weight: 148 pounds (67.1 kg)
- Shipping dimensions: 23"W x 29"D x 52"H (58.4 cm W x 73.7 cm D x 132.1 cm H)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



View similar products online at <https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300