



CAST CARTS

6025-TC

STAINLESS STEEL CAST CART

Eight Drawers with a Deluxe Accessory Package

Features:

- Stainless steel cabinet and drawers
- Painted aluminum top compartment
- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet/ pump pole
- 24" stainless steel ruler etched on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free
- Also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts.

Overall Cart Dimensions:

- 43" H x 44" W x 21" D (109.2 cm H x 111.8 cm W x 53.3 cm D) w/ top compartment (TC) 66" H x 44" W x 21" D (167.6 cm H x 111.7 cm W x 53.3 cm D)
- Cart net weight: 226 lbs (102.7 kg)
- Approximate shipping weight: 271 lbs (123.2 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D (144.8 cm H x 119.4 cm W x 68.6 cm D)



Engineered & Manufactured in the U.S.



Drawer Configuration:

- (1) - 3-inch with black vinyl mat (7.6 cm)
- (5) - 3-inch with 4 adj. dividers per drawer (7.6 cm)
- (2) - 6-inch with 4 adj. dividers per drawer (22.9 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Shown in stainless steel. More than 20 additional colors available in heat sealed powder coat steel.

Stainless Steel Eight Drawer Cast Cart with Deluxe Accessory Package



View similar products online at
<https://www.harloff.com/product-category/cast-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300