



# CAST CARTS

**6025E**

## STAINLESS STEEL CAST CART

Eight Drawers, Deluxe Accessory Package, E-Lock

### Features:

- Stainless steel cabinet and drawers
- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet/ pump pole
- 24" stainless steel ruler etched on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free
- Also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks (shown)
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts.



**Engineered & Manufactured in the U.S.**

**Shown in stainless steel. More than 20 additional colors available in heat sealed powder coat steel.**

#### Overall Cart Dimensions:

- 43" H x 38" W x 21" D (109.2 cm H x 96.5 cm W x 53.3 cm D) w/ pump pole
- 76" H x 44" W x 21" D (193 cm H x 111.8 cm W x 53.3 cm D)
- Writing surface: 42" H x (106.7 cm)
- Cart net weight: 207 lbs (94.1 kg)
- Approximate shipping weight: 252 lbs (114.5 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D (144.8 cm H x 119.4 cm W x 68.6 cm D)



### Drawer Configuration:

- (1) - 3-inch with black vinyl mat (7.6 cm)
- (5) - 3-inch with 4 adj. dividers per drawer (7.6 cm)
- (2) - 6-inch with 4 adj. dividers per drawer (22.9 cm)

### Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

**Stainless Steel Eight Drawer Cast Cart with Deluxe Accessory Package, E-Lock**



View similar products online at  
<https://www.harloff.com/product-category/cast-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)**  
**[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**