

Painted Steel Cast Cart/ Model: 6030

Eight Drawer, Standard Package



Over 20 Available Colors:

Beige, Burgundy, Light Blue, Light Gray, Navy, Red, Teal, Yellow, Cream, Taupe, White, Brushed Silver, Forest Green, Sand, Hammertone Black, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammertone Red, Purple, Mauve and Pink.

**Hammertone Red, Purple, Mauve & Pink are SPECIAL COLORS and if products are ordered in a special color they can not be returned.*

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Stainless steel bucket holder with vinyl mat
- 24" stainless steel ruler mounted on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free
- Designed in conjunction with orthopedic technicians to maximize functionality

Other Cast Cart Options:

- All cast carts are also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks. Call Harloff for details.

Drawer Configuration:

- (1) - 3-inch with black vinyl mat (7.6 cm)
- (5) - 3-inch with 4 adj. dividers per drawer (7.6 cm)
- (2) - 6-inch with 4 adj. dividers per drawer (22.9 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Overall Cart Dimensions:

- 43" H x 38" W x 21" D (109.2 cm H x 96.5 cm W x 53.3 cm D)
- Writing surface: 42" H x (106.7 cm)
- Cart net weight: 171 lbs (77.7 kg)
- Approximate shipping weight: 216 lbs (98.2 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D (144.8 cm H x 119.4 cm W x 68.6 cm D)