



# MEDICATION CARTS

## AL16CS-5

### CASSETTE MEDICATION CART 16 Bin Cassette Medication with a Key Lock

#### Features:

- Painted aluminum cabinet and drawers
- 16 each 5" wide (12.7cm) patient bin drawers in 2 removable cassettes; bin dividers and labels are included
- Cassettes are located in pull-out drawers for ergonomic removal
- BEST® lock key locking system with two keys
- Pull-out shelf on right side of cart
- Pre-threaded accessory mounting holes
- Self-closing, ball-bearing full extension drawer slides
- Removable plastic top with integrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one directional, two with brake
- Latex free
- Heat sealed powder coat finish



Color shown is Cream. More than 20 additional colors are available.

#### Suggested Cart Options:

- Cassette exchange package: 16CS5-PKG (includes 16 each 5" wide (12.7 cm) patient bin drawers and 2 double-tier aluminum cassettes)
- Narcotics drawer: CS-4FLN (separately locking narcotics drawer for 4" multi-purpose storage drawer)

#### OptimAL Line Cart Model: 16 Bin Cassette Medication with Key Lock, Cassette Exchange Package

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

#### Drawer Configuration:

- (1) - 4-inch (10.2 cm) multi-purpose storage drawer
- (2) - 8-inch (20.3 cm) double tier cassettes, with eight 5-inch wide patient bin drawers each
- (1) - 4-inch (20.3 cm) multi-purpose storage drawer
- (1) - 8-inch (20.3 cm) multi-purpose storage drawer
- Internal Drawer Dims: 24" W x 17" D (50.8 cm x 43.2 cm)
- Bin Dims: 5" W x 9.5" L x 3" H (12.7 cm x 24.1 cm x 7.6 cm)

#### Overall Cart Dimensions:

- 45.7" H x 30.5" W x 23.8" D (116.2 cm x 77.5 cm x 60.6 cm)
- Writing surface: 43.75" H (111.1 cm)
- Cart net weight: 132 lbs (59.9 kg)
- Approximate shipping weight: 180 lbs (81.8 kg)
- Approximate shipping dims: 48" H x 33" W x 26" D (121.9 cm x 83.8 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

**Engineered & Manufactured in the U.S.**



View similar products online at  
<https://www.harloff.com/product-category/med-bin-medication-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300