



MEDICATION CARTS

AL602PC

ALUMINUM PUNCH CARD MED CART 750 Punch Card Capacity with a Key Lock

Features:

- Dual column medication cart for punch card and liquid storage with a maximum capacity of 750 punch cards
- Includes one single width locking narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted aluminum cabinet and drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective gray corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Adjustable punch card row dividers



Color shown is Light Gray.
More than 20 additional colors
are available.

Aluminum 750 Punch Card Capacity Medication Cart with BEST® Lock on Cabinet, Standard Package

OptimAL Line medication carts for nursing homes are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.

Left Column Drawer Configuration:

- 1 ea - 4" multipurpose storage drawer
- 3 ea - 9.75" punch card drawers

Right Column Drawer Configuration:

- 1 ea - 4" multipurpose storage drawer
- 1 ea - 9.75" locking narcotics drawer
- 2 ea - 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm x 111.1 cm x 61 cm)
- Weight: 145 lbs (65.9 kg)
- Approx shipping weight: 190 lbs (86.4 kg)
- Approx Shipping Dims: 56" H x 46" W x 26" D (142.2 cm x 116.8 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered &
Manufactured
in the U.S.**



View similar products online at
<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**