



# MEDICATION CARTS AL8W10K4MOT

## MEDICINE-ON-TIME® MED CART 270 MOT® Card Capacity with a Key Lock

### Features:

- Single column medication delivery cart for MOT® cards and liquid storage with a maximum capacity of 270 Medicine-On-Time® cards
- Includes one single width locking narcotics box made for MOT® packaging that secures with an Illinois lock
- Key locking cart with 4 keys per lock
- Durable scratch resistant heat sealed powder coat finish
- All aluminum construction
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Removeable plastic top
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Painted steel top rail / push handle
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- Body has pre-threaded accessory mounting holes
- Latex free cart is delivered fully assembled



Color shown is Brushed Silver. More than 20 additional colors are available.

### Aluminum 270 Medicine-On-Time® Card Capacity Medication Cart with Key Lock on Cabinet

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

#### Internal Drawer Dimensions:

- 30" W x 17" D

#### Assembled Cart Dimensions:

- Writing Surface: 44"
- Cart weight: 125 pounds
- 48.5" H x 38.25" W x 25" D (123.2 cm x 97.2 cm x 63.5 cm)
- Approx shipping weight: 151 pounds
- Approx shipping dims: 56" H x 46" W x 26" D (142.2 cm x 116.8 cm x 66 cm)

#### Drawer Configuration:

- 1 ea – 3.25" multipurpose storage drawer
- 1 ea – 10" drawer with one aluminum riveted divider and locking narcotics box.
- 2 ea – 10" drawer with two aluminum riveted dividers.

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered & Manufactured in the U.S.**



View similar products online at  
<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300