



Slides Under Beds and Countertops

ETC-4

ETC LINE TREATMENT CART with Four Drawers and a Standard Key Lock

Harloff's smallest cart, this economy treatment cart is shorter than standard counter heights and is designed to slide under equipment or beds.

Features:

- Sturdy, uni-body steel construction
- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Latex free
- Heat sealed powder coat finish
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

Color shown is Cream.
More than 20 additional
colors are available.



Engineered &
Manufactured
in the U.S.



Four Drawer Treatment Cart with a Standard Key Lock

Drawer Configuration:

- (1) - 2.75-inch (6.98 cm)
- (3) - 5.75-inch (14.6 cm)

Assembled Cart Dimensions:

- Internal Drawer Dimensions: 20" W x 12.75" D (50.8 cm x 32.4 cm)
- 28.63" H x 23.5" W x 15" D (72.73 cm H x 59.4 W x 38.4 D)
- Cart weight: 55 lbs (25 kg)
- Approx. shipping weight: 86 pounds (39.1 kg)
- Approx. Shipping Dims: 33" H x 29" W x 23" D (83.8 cm H x 73.7 cm W x 58.4 cm D)

View similar products online at
<https://www.harloff.com/m-series-carts/>

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300