



M-series

M3DS1824E05

M-SERIES SHORT MINI CART

Five Drawers with an Electronic Pushbutton Lock

Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 24" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), one 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Color shown in Teal. More than 20 additional colors are available.

All dimensions are available in the "MA-Dims-18wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Engineered & Manufactured in the U.S.

Short Height, Mini Width, Five Drawer Anesthesia Cart with Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**