



Wood Style

MDS1830KC7-CM

M-SERIES TALL WOOD STYLE CART

with Seven Drawers and a Key Lock

Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.

Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Sturdy, uni-body steel construction
- Vinyl coated aluminum drawer fronts
- Choice of lock options:
www.harloff.com/locking-options/
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Removable plastic top
- Industry only fully curved protective bumper
- 30" of vertical drawer space
- Drawer Configuration: Five 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)
- Durable scratch resistant hard baked powder coat finish
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Pre-threaded mounting holes for accessories found at harloff.com/product-category/accessories/
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Engineered &
Manufactured
in the U.S.



All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

Colors:

Wood Style drawer fronts are shown in Cherry Mahogany, but are available in three standard options. Shown with a Light Gray body, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

Tall Height, Mini Width, Seven Drawer Wood Style Cart with a Standard Key Lock

Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of M-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff metal cart.



View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300