



M-series MDS3030B09PED

M-SERIES TALL PEDIATRIC CART with Nine Drawers & Breakaway Drawer Locks

Features:

- Sturdy, uni-body steel construction
- Tall cabinet with separate breakaway locks on each drawer
- Cart is designed to be compatible with popular pediatric color coding systems
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Ergonomic side push handles can be mounted on either side
- Comfort drawer pulls
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Painted steel breakaway lock bar
- Pack of 100 breakaway seals
- 30" of vertical drawer space
- Drawer Configuration: Eight 3" (7.6 cm) drawers and one 6" (15.2 cm) drawer
- Latex free
- Cart is delivered assembled
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



**Color shown is Beige with color coded drawers.
More than 20 additional body colors are available.**

Tall, Standard Width, Nine Drawer Pediatric Cart

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



**Engineered &
Manufactured
in the U.S.**



View similar products online at
<https://www.harloff.com/product-category/pediatric-resuscitation-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**