




M-series MDS3030B09PED

M-SERIES TALL PEDIATRIC CART with Nine Drawers & Breakaway Drawer Locks

Features:

- Sturdy, uni-body steel construction
- Tall cabinet with separate breakaway locks on each drawer
- Cart is designed to be compatible with popular pediatric color coding systems
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Ergonomic side push handles can be mounted on either side
- Comfort drawer pulls
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Painted steel breakaway lock bar
- Pack of 100 breakaway seals
- 30" of vertical drawer space
- Drawer Configuration: Eight 3" (7.6 cm) drawers and one 6" (15.2 cm) drawer
- Latex free
- Cart is delivered assembled
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

Color shown is Beige with color coded drawers. More than 20 additional body colors are available.



**Engineered &
Manufactured
in the U.S.**

Tall, Standard Width, Nine Drawer Pediatric Cart

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at
<https://www.harloff.com/product-category/pediatric-resuscitation-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**