



# MDS3030K06+MD30-ADJSHLF2-XBAR

with Adjustable Overhead Shelf

## M-SERIES TALL MEDICAL EQUIPMENT CART

with MD30-ADJSHLF2-XBAR Accessory

### Tall, Standard Width, Six Drawer Monitor Cart with Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

### Accessories Include:

Adjustable Overhead Shelf (MD30-ADJSHLF2-XBAR)



### Features:

- Sturdy, uni-body steel construction
- Choice of lock options: [www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Removable plastic top
- Fully customizable drawer configurations
- Ergonomic side push handles can be mounted on either side
- 30" of vertical drawer space
- Cart is delivered assembled
- Latex free

## M-series

Color shown is Navy. More than 20 additional colors are available.

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300  
[medicalsales@harloff.com](mailto:medicalsales@harloff.com) • [www.harloff.com](http://www.harloff.com)