



# MDS3030KC4+MD30-DRW9-EH

## Equipment Holder

**M-SERIES STANDARD WIDTH CART**  
with 9" Equipment Holder, 3 Drawers and Key Lock

### M-Series Features:

- Sturdy, uni-body steel construction
- Choice of lock options:  
[www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)
- Up to 30" of fully customizable vertical space is available for both drawer and equipment holder configurations
- 17" deep open equipment holding area can replace standard drawer areas in three possible height options
- 2" pass-through hole in the rear for cords and cables
- Full extension drawer slides for easy access
- Cart is delivered assembled
- Removable plastic top
- Industry only fully curved protective bumper
- Durable scratch resistant hard baked powder coat finish
- Body has pre-threaded accessory mounting holes
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free

### Colors:

Shown in Red, but available in more than twenty color options which you can view online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/).



## M-series

#### Drawer Configuration:

- One MD30-DRW9-EH equipment holder with an 8" H x 24.625" W x 17" D (20.3 cm H x 62.5 cm W x 43.2cm D) opening, one 3" drawer (7.6 cm) and two 9" drawers (22.9 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "[MA-Dims-30wide-5caster.pdf](#)" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

Designed and built to meet the needs of today's healthcare market. Replace 9", 12 or 18" of vertical drawer space at the top or bottom of your M-Series cart with a doored equipment holder, perfect for storing electronics or large medical equipment. M-Series carts offer you choices.

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300