



# MALIGNANT HYPERTHERMIA

## MH4200B

### MALIGNANT HYPERTHERMIA CART with 1.8 Cubic Feet Medical Grade Refrigerator

#### Features:

- 20-gauge powder-coated steel construction w/built-in Follett® REF2 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- Highest quality in the industry Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Breakaway locks on drawers and refrigerator (refrigerator can be keylocked)
- Cart ships ready to use with wheels installed
- Removable plastic top

#### Colors:

Shown in Hammertone Blue, but available in over twenty color options. View them online: [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/)



### Malignant Hyperthermia Cart, Two Drawer, Breakaway Lock with No Accessories

#### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 48.5" H x 37.5" W x 26" D (123.2 cm H x 95.3cm W x 66cm D)
- Cart net weight: 188 lbs (85.3 kg)
- Approximate shipping weight: 227 lbs (103.2 kg)
- Approximate shipping dims: 62" H x 44" W x 30" D (157.5cm H x 111.8cm W x 76.2cm D)

#### Drawer Configuration:

- (2) 3" (7.6 cm) and (1) Compartment with pre-installed refrigerator



Engineered &  
Manufactured  
in the U.S.

#### Optional Accessories:



Large wire basket



Side storage  
cabinets



Two large chart/  
document holders



Waste container &  
mounting bracket



View similar products online at  
<https://www.harloff.com/malignant-hyperthermia-cart/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300