




**MPA1821K03**

**ALUMINUM**

## **A-SERIES MINI WIDTH CART** with Three Drawers and a Standard Key Lock

### **Features:**

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters
  - one steering caster and two braking casters
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <https://www.harloff.com/locking-options/>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 21" of vertical drawer space
- Drawer Configuration: Two 6" drawers (15.2 cm) and one 9" drawer (22.9 cm)
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Color shown with a Taupe body and Cream drawers. More than 20 additional colors are available.



**Engineered & Manufactured in the U.S.**

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

### **X-Short Height, Mini Width, Three Drawer Cart with a Standard Key Lock**

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Light-weight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



View similar products online at  
<https://www.harloff.com/product-category/a-series-medical-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)**  
**[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**