



series

ALUMINUM

A-SERIES MINI WIDTH CART

with Five Drawers and a Breakaway Lock

Features:

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
 - Curved drawer pulls, durable and elegant
 - Fully welded uni-body aluminum construction, highest quality in the industry
 - Heat-sealed durable scratch resistant powder coat finish, inside and out
 - Ordered accessories are easily attached to pre-threaded holes with a screwdriver
 - Industry only fully curved protective bumper guards your cart and walls from damage
 - Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
 - Removable plastic top
 - Fully customizable drawer configurations (no extra charge) with 24" of vertical drawer space
 - Industry best 5" quiet rolling casters - one steering caster and two braking casters
 - Multiple alternative storage options for items that won't fit inside a traditional drawer
 - Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.9 cm)
 - Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems.
- View online at <https://www.harloff.com/locking-options/>
- Cart is delivered assembled
 - 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Engineered & Manufactured in the U.S.



Color shown with a Red body and White drawers. More than 20 additional colors are available.

Optional MD18-EMG Accessory Package (shown) Includes:



Oxygen tank holder (O2HLDR)



Cardiac board and brackets (MD-CARDBRD)



Aluminum 16" Push Side Rail (SIDERAIL16)



IV pole (MD-IV2)



All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300