



MPA3030B09PED

ALUMINUM

A-SERIES TALL PEDIATRIC CART **with Nine Drawers and Breakaway Drawer Locks**

Features:

- Designed for treating pediatric emergencies by using a multi-colored drawer labeling system comparable to the Broselow color coding method
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top, easily replaceable
- Ergonomic side push handles can be mounted on the left or right
- Painted steel break-away lock bar on every drawer
- Pack of 100 break-away seals
- Industry best 5" quiet rolling casters
 - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space
- Drawer Configuration: Eight 3" (7.6 cm) drawers and one 6" (15.2 cm) drawer
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Color shown with a Beige body and color coded drawers. More than 20 additional body colors are available.

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Engineered & Manufactured in the U.S.

Tall, Standard Width, Nine Drawer Pediatric Cart with Breakaway Drawer Locks

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



View similar products online at

<https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300