



MPA3030WLP16

ALUMINUM

A-SERIES STANDARD WIDTH CART

Six Drawer, Wireless Elec. Lock w/Keypad & Prox

Features:

- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top, easily replaceable
- Ergonomic side push handles can be mounted on the left or right
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters
 - 1 steering caster and 2 braking casters
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <https://www.harloff.com/locking-options/>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space
- Drawer Configuration: Four 3" drawers (7.6 cm), one 6" drawers (15.2 cm) and one 12" drawer (30.4 cm)
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



**Color shown
in Light Blue.
More than
20 additional
colors are
available.**

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



**Engineered &
Manufactured
in the U.S.**

Tall, Standard Width, Six Drawer Cart and Wireless Electronic Lock with Keypad and Proximity Reader

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



View similar products online at

<https://www.harloff.com/product-category/a-series-medical-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**