



# MEDICATION CARTS **SL100DSPN**

## **DISPILL® MEDICATION CART** 144 Dispill® Card Capacity with a Key Lock

### Features:

- Single column medication delivery cart for Dispill® cards and liquid storage with a maximum capacity of 144 Dispill® cards
- Key lock system with Best® brand lock on cabinet
- Illinois lock on narcotics box located in top 9" drawer, designed to hold Dispill® cards
- Painted steel cabinet with aluminum drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective gray corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Adjustable punch card row dividers
- Pre-threaded accessory mounting holes



Color shown is Light Blue. More than 20 additional colors are available.

### 144 Dispill® Card Capacity Medication Cart with BEST® Lock on Cabinet

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

#### Drawer Configuration

- 1 ea - 4" multipurpose storage drawer
- 1 ea - 3.25" multipurpose storage drawer
- 1 ea - 13" Dispill card drawer with dividers and locking narcotics box
- 1 ea - 13" Dispill card drawers with dividers

#### Assembled Cart Dimensions:

- Writing surface: 44" (111.8 cm)
- Cart weight: 175 lbs (79.4 kg)
- 45.75" H x 36" W x 23.5" D (116.2 cm x 91.4 cm x 59.7 cm)
- Approx shipping weight: 215 lbs (97.5 kg)
- Approx shipping dims: 57" H x 39" W x 26" D (144.8 cm x 99.1 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered & Manufactured in the U.S.**



View similar products online at  
<https://www.harloff.com/product-category/specialty-packaging-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)**  
**[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**