



MEDICATION CARTS

SL28BIN5

STANDARD LINE MED BIN CART with 28 – 5 inch Bins, Key Locking

Features:

- Painted steel cabinet and aluminum drawers
- Dual column medication cart with 28 - 5" wide medication bins
- Dividers and labels included
- Single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Heat sealed powder coat finish
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts



Color shown is Hammertone Blue. More than 20 additional colors are available.

Standard Line Med Bin Cart with 28 - 5" Bins, Key Locking

Harloff medication carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care.

Left Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 7 ea – 3.25" med bin drawers with 4 – 5" bins per drawer
- 1 ea – 6.5" multipurpose drawer with adjustable dividers

Right Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.625" H x 43.5" W x 24" D (115.9 cm H x 110.5 cm W x 61 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 235 pounds (106.8 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 68.6 cm D)
- Approx. shipping weight: 302 pounds (137 kg)



Engineered &
Manufactured
in the U.S.

View similar products online at

<https://www.harloff.com/product-category/specialty-packaging-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300

HARLOFF
INNOVATING STORAGE FOR CARE