



# MEDICATION CARTS

## SL294BOX

### UNIT-DOSE MEDICATION BOX CART

294 Box Capacity, 2-Column, Standard Line

#### Features:

- Painted steel cabinet and aluminum drawers
- Dual column medication cart with 294 medication box capacity, single width narcotics drawer
- BEST® brand lock on cabinet with four keys
- Illinois lock on narcotics drawer with four keys
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Pull-out side shelf
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Heat sealed powder coat finish
- Dividers are available to accommodate different medication packaging
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Color shown is Brushed Silver. More than 20 additional colors are available.



#### Left Column Drawer Configuration

- 1 ea - 4" multipurpose storage drawer
- 6 ea - 3.25" med box drawers with internal dividers, capacity 42 boxes per drawer
- 1 ea - 9.75" multipurpose storage drawer with adjustable dividers

#### Right Column Drawer Configuration

- 1 ea - 4" multipurpose storage drawer
- 1 ea - 9.75" locking narcotics drawer
- 2 ea - 9.75" multipurpose storage drawer

#### Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 210 pounds (95.5 kg)
- Approx. shipping weight: 255 pounds (116 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm x 119.4 cm x 68.6 cm)



Engineered & Manufactured in the U.S.

## Standard Line 294 Capacity Unit Dose Medication Box Cart

Harloff med carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care.



View similar products online at

<https://www.harloff.com/product-category/multi-dose-medication-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)

[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300