



MEDICATION CARTS

SL36BIN5

STANDARD LINE MED BIN CART with 36 – 5 inch Bins, Key Locking

Features:

- Painted steel cabinet and aluminum drawers
- Dual column medication cart with 36 - 5" wide medication bins
- Dividers and labels included
- Single width narcotics drawer
- CompX® electronic lock on cabinet with BEST® brand key lock override
- Illinois lock on narcotics drawer
- Access right may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set of programmed time)
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Heat sealed powder coat finish
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts



Standard Line Med Bin Cart with 36 - 5" Bins, Key Locking

Harloff medication carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care.

Left Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 9 ea – 3.25" med bin drawers with 4 – 5" bins per drawer

Right Column Drawer Configuration:

- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 215 pounds (97.7 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 68.6 cm D)
- Approx. shipping weight: 270 pounds (122.7 kg)



**Engineered &
Manufactured
in the U.S.**



View similar products online at

<https://www.harloff.com/product-category/specialty-packaging-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**