



# MEDICATION CARTS

## SL550PC

### PUNCH CARD MEDICATION CART 600 Punch Card Capacity with a Key Lock

#### Features:

- Single column medication cart for punch card and liquid storage with a maximum capacity of 600 punch cards
- Includes one in-drawer single width locking narcotics box
- BEST® brand lock (with four keys) on cabinet; Illinois lock (with four keys) on narcotics drawer
- Painted steel cabinet with aluminum drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Protective gray corner bumpers
- Plastic replaceable top with integrated push handles
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Adjustable punch card row dividers



Engineered & Manufactured in the U.S.

Color shown is Light Gray. More than 20 additional colors are available.



#### Standard 600 Punch Card Capacity Medication Cart with BEST® Lock on Cabinet, Standard Package

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

##### Left Column Drawer Configuration:

- (1) – 4" multipurpose storage drawer
- (1) – 9.75" punch card drawer with dividers and locking narcotics box
- (2) – 9.75" punch card drawers with dividers

##### Overall Cart Dimensions:

- Writing surface: 44" (111.8 cm)
- 45.75" H x 36" W x 23.5" D (116.2 cm x 91.4 cm x 59.7 cm)
- Cart weight: 170 lbs (77.3 kg)
- Approx shipping weight: 220 lbs (100 kg)
- Approx shipping dims: 57" H x 36" W x 26" D (144.8 cm x 99.1 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at

<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300