



MED CARTS

SL600PC-EKC

PUNCH CARD MEDICATION CART 750 Punch Card Capacity w/Elec. Keypad Lock

Features:

- Dual column medication cart for punch card and liquid storage with a maximum capacity of 750 punch cards
- Electronic locking cart with keypad access and key lock override with a BEST® brand lock
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart. May add, edit, or delete user and user groups
- Operates with six D-cell Batteries; life 1 year minimum at 25-50 uses/day
- Lock has an LED low battery indicator
- Includes one single width locking narcotics drawer secured by an Illinois lock
- Durable scratch resistant heat sealed powder coat finish steel cabinet with aluminum drawers
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective gray corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Adjustable punch card row dividers



Standard 750 Punch Card Capacity Medication Cart with Electronic Keypad Lock, Standard Package

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

Left Column Drawer Configuration:

- 1 ea - 4" multipurpose storage drawer
- 3 ea - 9.75" punch card drawers

Right Column Drawer Configuration:

- 1 ea - 4" multipurpose storage drawer
- 1 ea - 9.75" locking narcotics drawer
- 2 ea - 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- Cart weight: 195 lbs (88.6 kg)
- 45.75" H x 43" W x 23.5" D (116.2 cm x 109.2 cm x 59.7 cm)
- Approx shipping weight: 242 lbs (106.8 kg)
- Approx Shipping Dims: 58" H x 46" W x 26" D (147.3 cm x 116.8 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Engineered & Manufactured in the U.S.



View similar products online at <https://www.harloff.com/product-category/punch-card-carts-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300