



# MEDICATION CARTS

## SL72MD

### MULTI-DOSE BIN MEDICATION CART 72 Bin Capacity, 2-Column, Standard Package

#### Features:

- Painted steel cabinet and aluminum drawers
- Dual column medication cart with 72 multi-dose bin capacity
- Dividers and labels included
- Single width narcotics drawer with Illinois lock
- BEST® brand lock on cabinet
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Heat sealed powder coat finish
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Color shown is Brushed Silver. More than 20 additional colors are available.



#### Standard Line Multi-Dose Bin Medication Cart, Standard Package

Harloff med carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care.

##### Left Column Drawer Configuration

- 6 ea - 4" drawers with internal dividers, capacity 12 bins per drawer (Bins are not included)
- 1 ea - 9.75" multipurpose storage drawer w/adjustable dividers

##### Right Column Drawer Configuration

- 2 ea - 4" multipurpose storage drawer
- 1 ea - 8" locking narcotics drawer
- 2 ea - 8" multipurpose storage drawer

##### Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 220 pounds (100 kg)
- Approx. shipping weight: 265 pounds (120.5 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm x 119.4 cm x 68.6 cm)



Engineered & Manufactured in the U.S.



View similar products online at

<https://www.harloff.com/product-category/multi-dose-medication-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300