



MEDICATION CARTS VLT150EMOT

MEDICINE-ON-TIME® MED CART 150 MOT® Card Capacity with an E-Lock

Features:

- Single column medication delivery cart for MOT® cards and liquid storage with a maximum capacity of 150 Medicine-On-Time® cards
- Includes one single width locking narcotics box made for MOT® packaging that secures with an Illinois lock
- Basic electronic pushbutton lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Painted steel cabinet with aluminum drawers
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Painted steel top rail / push handle
- Removeable plastic top
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Protective gray corner bumpers
- Body has pre-threaded accessory mounting holes
- Latex free cart is delivered fully assembled



Color shown is Light Gray. More than 20 additional colors are available.

Value 150 Medicine-On-Time® Card Capacity Medication Cart with Electronic Pushbutton Lock on Cabinet

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawers
- 1 - ea 9.75" drawer with locking narcotics compartment
- 2 - ea 9.75" drawers with one adjustable divider

Assembled Cart Dimensions:

- Cart weight: 140 pounds (63.7 kg)
- 46.75" H x 29.5" W x 22.5" D (121.3 cm H x 75 cm W x 57.2 cm D)
- Approx. shipping weight: 180 pounds (81.8kg)
- Approx Shipping dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.



View similar products online at
<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300