



# MEDICATION CARTS VLT216EBOX

## VALUE LINE UNIT-DOSE MED CART with Narc Box and Electronic Pushbutton Lock

### Features:

- Painted steel cabinet and drawers
- Tall cabinet with basic electronic pushbutton lock, key lock override with four keys
- Factory set user and supervisor code; lock can be programmed in field
- Internal locking narcotics box in bottom 9" drawer, 4 keys
- Capacity of up to 216 unit dose boxes or 60 sections for strip packaging
- Pull out Side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- Painted steel top rail/push handle
- Plastic top, removable
- Full extension ball-bearing drawer slides
- 5" ball-bearing premium casters, all full swivel – one with brake, one directional
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Heat sealed powder coat finish
- Certified "Latex Free"
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts



Engineered &  
Manufactured  
in the U.S.



**Color shown is Hammertone Gray with Light Gray Drawers. More than 20 additional colors are available.**

#### Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

#### Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 6 - ea 3" drawers with internal unit dose dividers
- 1 - ea 9" drawer with 1 divider and locking narcotics box

#### Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 170 pounds (77.3 kg)
- Approx. shipping weight: 210 pounds (95.5kg)

### Value Line Unit-Dose Med Cart with Narc Box & E-Lock

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



View similar products online at

<https://www.harloff.com/product-category/med-bin-medication-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**