



# MEDICATION CARTS

# VLT480EPC

## PUNCH CARD MEDICATION CART 600 Punch Card Capacity w/Pushbutton E-Lock

### Features:

- Single column medication cart for punch card and liquid storage with a maximum capacity of 600 punch cards
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Includes one single width locking narcotics drawer with Illinois lock
- Painted steel cabinet with aluminum drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Painted steel top rail / push handle
- Removeable plastic top
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Protective gray corner bumpers
- Adjustable punch card row dividers
- Body has pre-threaded accessory mounting holes
- Cart is delivered fully assembled
- Latex free



Color shown is White. More than 20 additional colors are available.

### Value 600 Punch Card Capacity Med Cart with Electronic Pushbutton Lock on Cabinet, Standard Package

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

#### Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 1 - ea 9.75" drawer with one divider and locking narcotics box
- 2 - ea 9.75" drawers with two adjustable dividers

#### Assembled Cart Dimensions:

- Cart weight: 155 pounds (70.45 kg)
- 46.75" H x 33" W x 22.5" D (121.3 cm H x 82.5 cm W x 57.2 cm D)
- Approx. shipping weight: 185 pounds (89 kg)
- Approx Shipping dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered &  
Manufactured  
in the U.S.**



View similar products online at  
<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**