



Wood Vinyl

WV240PC-CM

PUNCH CARD MEDICATION CART

240 Capacity Wood Vinyl Covered Aluminum

Features:

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Single column medication cart with 240 punch card capacity
- Key locking cabinet with locking narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- 5" ball-bearing premium rubber casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Drawer Configuration:

- 1 ea – 3.1" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics box and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers

Skirt Option:

- #WV-2WSKIRT-CM – Skirt option to cover wheels

Colors:

Shown in Cherry Mahogany. Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



(CM) Cherry Mahogany



Overall Cart Dimensions:

- 44.5" H x 19" W x 25" D (113 cm x 48.3 cm x 63.5 cm)
- Writing surface: 44" (111.8 cm)
- Cart weight: 102 pounds (46.3 kg)
- Approx. shipping weight: 135 pounds (61.2 kg)

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.



Engineered & Manufactured in the U.S.



View similar products online at

<https://www.harloff.com/product-category/medication-furniture-carts-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com

www.harloff.com • Ph: (719) 637-0300