



Wood Vinyl

WV360PC-CM

PUNCH CARD MEDICATION CART

360 Capacity Wood Vinyl Covered Aluminum

Features:


- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Single column medication cart with 360 punch card capacity
- Key locking cabinet and narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 5" ball-bearing premium rubber casters, full swivel, one with toe locking brake and one directional
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Drawer Configuration:

- 1 ea – 3" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics box and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers

Skirt Option:

- #WV-3WSKIRT-CM – Skirt option to cover wheels



(CM) Cherry Mahogany

Colors:

Shown in Cherry Mahogany color. Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



Engineered & Manufactured in the U.S.

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.

Overall Cart Dimensions:

- 44.5" H x 30" W x 23.5" D (113.1 cm H x 76.2 cm W x 59.7 cm D)
- Writing Surface: 44" (111.8)
- Cart weight: 112 lbs (50.8 kg)
- Approximate shipping weight: 142 lbs (64.4 kg)

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/medication-furniture-carts-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300