



*Wood Vinyl*

**WV540PC-CM**

## **PUNCH CARD MEDICATION CART**

540 Capacity Wood Vinyl Covered Aluminum

### **Features:**

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Single column medication cart with 540 punch card capacity
- Key locking cabinet and narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

### **Drawer Configuration:**

- 1 ea – 3" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics box and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers

### **Skirt Option:**

- #WV-3WSKIRT-CM – Skirt option to cover wheels

### **Colors:**

Shown in Cherry Mahogany color. Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



(CM) Cherry Mahogany

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.

#### **Overall Cart Dimensions:**

- 45" H x 34.5" W x 24.5" D
- Writing Surface: 44" (111.8)
- Cart weight: 127 lbs (57.7 kg)
- Approximate shipping weight: 157 lbs (71.2 kg)



**Engineered & Manufactured in the U.S.**

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/medication-furniture-carts-cabinets/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)

[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300