



*Wood Vinyl*

**WV6400SC-CM**

## TALL ALUMINUM EMERGENCY CART

Wood Vinyl with Side Cabinet and Breakaway Lock

### Features:


- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Breakaway lock bar (wood look) with 100 breakaway seals
- Includes side cabinet (wood look inside and out) with shelving
- Pull out shelf
- Ball-bearing, full extension drawer slides
- 5" (12.7cm) ball-bearing premium casters, all full swivel, one braking, one directional
- This cart is designed to look like a piece of furniture
- Cart is washable with soap and water
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

### Skirt Option:

- #WV-4WSKIRT-CM – Skirt option to cover wheels

### Colors:

Shown in Cherry Mahogany. Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



(CM) Cherry Mahogany



**Engineered & Manufactured in the U.S.**

#### Drawer Configuration:

- Four 3.1" drawers (7.9 cm), one 6.4" drawer (16.4 cm) and one 13.1" drawer (33.3 cm)

#### Overall Cart Dimensions:

- 44.5" H x 39" W x 22.5" D (113 cm H x 99.1 cm W x 57.2 cm D)
- Writing Surface: 44" height (111.8 cm)
- Cart weight: 130 pounds (59 kg)
- Approx. shipping weight: 170 pounds (77.1 kg)

Bring the warm appearance of wood to your treatment and procedure carts with vinyl-coated aluminum carts. Designed and built to support culture change initiatives, this cart gives the buyer the best of both worlds, combining the aesthetics of wood with the long-term functionality of a Harloff metal cart.

View similar products online at

<https://www.harloff.com/product-category/medical-furniture-carts/>



**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
www.harloff.com • Ph: (719) 637-0300**