



Wood Vinyl

WV6520-CM

TALL 4 DRAWER TREATMENT CART

Wood Vinyl Coated Aluminum with Key Lock

Features:

- Lightweight vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate
- Key lock with two keys
- Pull out shelf
- Ball-bearing, full extension drawer slides
- 5" (12.7cm) ball-bearing premium casters, all full swivel, one braking, one directional
- This cart is designed to look like a piece of furniture
- Cart is washable with soap and water
- Twelve Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Skirt Option:


- #WV-3WSKIRT-CM – Skirt option to cover wheels



Engineered & Manufactured in the U.S.

Colors:

Shown in Cherry Mahogany. Harloff offers two standard vinyl color options, plus additional special order vinyl color options with a minimum cabinet purchase and extended lead time.



(CM) Cherry Mahogany



Drawer Configuration:

- One 3.1" drawer (7.9 cm) and three 9.75" drawers (24.8 cm)

Overall Cart Dimensions:

- 44.5" H x 30" W x 23.5" D (113 cm H x 76.2 cm W x 59.7 cm D)
- Writing Surface: 44" height (111.8 cm)
- Cart weight: 112 pounds (50.8 kg)
- Approx. shipping weight: 142 pounds (64.4 kg)

Bring the warm appearance of wood to your treatment and procedure carts with vinyl-coated aluminum carts. Designed and built to support culture change initiatives, this cart gives the buyer the best of both worlds, combining the aesthetics of wood with the long-term functionality of a Harloff metal cart.

View similar products online at

<https://www.harloff.com/product-category/medical-furniture-carts/>

HARLOFF
INNOVATING STORAGE FOR CARE

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300