



## For IV and Epidural Procedures

3123

### START CART WITH ARTICULATING ARM with Specialty Package, Four Drawers & Latch Lock

*Harloff's Start Cart is designed to reduce the risk of accidental needle sticks to patients and staff while improving the efficiency of the IV or epidural procedure.*



Engineered &  
Manufactured  
in the U.S.

#### Features:

- Heat sealed powder coat steel cabinet & drawers
- Pontoon bumper protects facility walls and cart
- Articulating arm allows clinician to work from sitting or standing position
- Dual lipped working surfaces with clear plastic cover overlay
- Adjustable push handle in 1-inch increments
- Soft grip handles for moving cart from sitting position
- 3 gallon sharps and waste containers make it easy to dispose of supplies
- Ballasted for stability and lower center of gravity
- 16" full extension drawer slides
- Key lock with two keys
- Premium 5-inch casters make it easy to roll around the OR and patient rooms
- Removable plastic top

#### Benefits:

- Clinicians will no longer have to gather needed equipment from a distant cart, resulting in going back and forth for supplies from the cart to the patient site.
- Clinicians will no longer have to use the patient's bed, stretcher or surgical operating room table as a work surface.
- Clinicians will no longer have to leave exposed needles on top of the patient's bed that may cause injury to the doctor, nurse, medical technician, or housekeeper
- Technicians can prepare supplies to be used in advance by clinicians which will save time and money.

**Color shown is Sand. More than 20 additional colors are available.**

Dual work surfaces and an articulating arm provide work area to prepare and dispose of IV and Epidural supplies. Cart provides compact storage for point of use supplies.



#### Assembled Cart Dimensions:

- Base cart: 21.50" w x 22.50" d (54.61 cm W x 57.15 cm D)
- Cart height Base Unit: 30" h (76.20 cm H)
- Articulating arm extends height to 45.75" h (116.21 cm H)
- Cart with waste and sharps containers: 31.50" w x 22.50" d (80.01 cm W x 57.15 cm D)
- Adjustable Push Handle: Lowest height 34.75" extends to 39.75" in 1" increments (Metric equivalent: Lowest height 88.27 cm extends to 100.97 cm in 2.54 cm increments)
- Cart net weight: 116 pounds (52.62 kg)
- Approximate shipping weight: 148 pounds (67.1 kg)
- Shipping dimensions: 23"W x 29"D x 52"H (58.4 cm W x 73.7 cm D x 132.1 cm H)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada.  
1-Year Limited Warranty on any bought-out parts.



View similar products online at

<https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)**  
**[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**