

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

Mini Line Anesthesia Cart/ Model: 3145K-ANS

Five Drawer with Key Lock, Specialty Accessory Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 24" (76.2 cm) of vertical drawer space
- Key lock with two keys
- 3" (7.6 cm) single wheel stem caster set, full swivel two with brake
- Pre-threaded accessory mounting holes
- · Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- · Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-tier raised back rail system (MINI-RBR)
- (2) Aluminum side rails/push handles (68161)
- 6-bin tilt bin organizer (683805)

Drawer Configuration:

- (3) 3-inch (7.6 cm)
- (1) 6-inch (15.2 cm)
- (1) 9-inch (22.9 cm)

Internal Drawer Dimensions:

• 12.6" W x 16.75" D (32 cm x 42.5cm)

Overall Cart Dimensions:

- 33.75" H x 19.25" W x 19.75" D (87.4 cm H x 49.5 cm W x 50.2 cm D)
- With accessories: 54.75" H x 21.25" W x 19" D
- Writing surface: 33.5" H x (81.2 cm)
- Cart net weight: 80 lbs (36.4 kg)
- Approximate shipping weight: 95 lbs (43.2 kg)
- Approximate shipping dims: 45" H x 24" W x 24" D (114.3 cm H x 60.9 cm W x 60.9 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts