

Mini Line Anesthesia Cart/ Model: 3145K

Five Drawer with Key Lock, Standard Package



Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5cm)

Overall Cart Dimensions:

- 34.4" H x 19.25" W x 19.75" D (87.4 cm H x 49.5 cm W x 50.2 cm D)
- Writing surface: 33.5" H x (81.2 cm)
- Cart net weight: 74 lbs (33.6 kg)
- Approximate shipping weight: 108 lbs (48.9 kg)
- Approximate shipping dims: 45" H x 24" W x 24" D (114.3 cm H x 60.9 cm W x 60.9 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts