

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

# Mini Line Emergency Cart/ Model: 3156B-EMG

Six Drawer with Breakaway Lock, Specialty Accessory Package



### Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### **Construction:**

Painted steel cabinet and drawers

#### Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 30" (76.2 cm) of vertical drawer space
- Painted steel break-away lock bar with 100 breakaway seals
- 3" (7.6 cm) single wheel stem caster set, full swivel two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

## Accessory Package Includes:

- 2-prong extendable IV pole (680401) ( IV pole extends to 76.5")
- Oxygen tank holder (AL680403)
- Cardiac board and brackets (MINI-CARDBRD)
- Aluminum push handle / side rail (68161)

# **Drawer Configuration:**

- (3) 3" (7.6 cm)(2) 6" (15.2 cm)
- (1) 9" (22.9 cm)

### **Internal Drawer Dimensions:**

12.6" W x 16.75" D (32 cm x 42.5 cm)

### **Overall Cart Dimensions:**

- 69.5" H x 24.5" W x 19.75" D (176.5 cm H x 62.2 cm W x 50.2 cm D) Height includes IV Pole – extends to 84.75" H (215.3 cm H)
- Writing surface: 39.75" H x (101 cm)
- Cart net weight: 107 lbs (48.6 kg)
- Approximate shipping weight: 136 lbs (61.8 kg)
- Approximate shipping dims: 52" H x 24" W x 24" D (106.7 cm H x 60.9 cm W x 60.9 cm D)

## Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts