



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Three Drawer Mini Line Infection Control Cart with Extended Stabilizing Bumper and 5" Wheels, with Key Lock #3243K

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints



Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 6-inch (15.2 cm)
- (2) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

Assembled Cart Dimensions:

- 37" H x 19.5" W x 25" D
(94.0 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 94 pounds (42.6 kg)
- Approximate shipping weight: 114 pounds (51.7 kg)

Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>