

Mini Line Emergency Cart/ Model: 3256B-EMG

Six Drawer with Breakaway Lock, Specialty Accessory Package

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics, and doctor offices with both space and budget constraints



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 30" (76.2 cm) of vertical drawer space
- Painted steel break-away lock bar with 100 break-away seals
- Stabilizing pontoon bumpers
- 5" (12.7 cm) ball-bearing caster set, all full swivel, two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-prong extendable IV pole (680401) (IV pole extends to 76.5")
- Oxygen tank holder (AL680403)
- Cardiac board and brackets (MINI-CARDBRD)
- Aluminum push handle / side rail (68161)

Drawer Configuration:

- (3) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 9" (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5 cm)

Overall Cart Dimensions:

- 68.25" H x 25" W x 22.2" D (173.4 cm H x 63.5 cm W x 56.4 cm D) Height includes IV pole in low position – extends to 86.25" (219.1 cm H)
- Writing surface: 41.75" H x (106 cm)
- Cart net weight: 103 lbs (46.8 kg)
- Approximate shipping weight: 130 lbs (59.1 kg)
- Approximate shipping dims: 52" H x 29" W x 24" D (131.1 cm H x 73.7 cm W x 60.9 cm D)