

Anesthesia Cart/ Model: 3256E-ANS

Six Drawer Mini Line Anesthesia Cart with Extended Stabilizing Bumper 5" Wheels, with Basic Electronic Pushbutton Lock, Specialty Accessory Package

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Stabilizing pontoon bumpers
- 5" (12.7 cm) ball-bearing caster set, all full swivel, two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-tier raised back rail system (MINI-RBR)
- (2) Aluminum side rails/push handles (68161)
- 6-bin tilt bin organizer (683805)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

Overall Cart Dimensions:

- 58.75" H x 19.5" W x 25" D (149.2 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 107 pounds (48.5 kg)
- Approximate shipping weight: 127 pounds (57.6 kg)
- Approximate shipping dims: 52" H x 29" W x 24" D (132 cm H x 73.6cm W x 60.9 cm D)