

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Three Drawer Mini Line Emergency Cart with Breakaway Lock, Standard Package #4143B

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices.



Color Shown is Red Additional Colors Available.

Shipping Dimensions:

- 49" H x 33" W x 26" D (124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 121 pounds (54.88 kg)

Specifications

Construction:

Painted steel cabinet and drawers

Features:

- Break-away locking system
- Space-saver design allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- 3" (7.6 cm) single-wheel stem caster set, full swivel two with brakes
- · Ball-bearing full extension drawer slides
- · Removable plastic top
- Pre-threaded accessory mounting holes
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

• 18.5" W x 16.75" D (47cm x 42.5cm)

Assembled cart dimensions:

- 33.75" H x 24" W x 18"D (85.7 cm H x 61 cm W x 45.7 cm D)
- Cart net weight: 75 pounds (34.1 kg)
- Approximate shipping weight: 82 pounds (37.19 kg)

Drawer Configuration:

- (1) 6-inch (15.2 cm)
- (2) 9-inch (22.9 cm)