

#### Stainless Steel Cast Cart/ Model: 6025

Eight Drawer, Deluxe Package



# Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

# **Construction:**

Stainless steel cabinet and drawers

#### Features:

- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet/ pump pole
- 24" stainless steel ruler etched on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free

# **Other Cast Cart Options:**

• All cast carts are also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks. Call Harloff for details.

#### **Drawer Configuration:**

- (1) 3-inch with black vinyl mat (7.6 cm)
- (5) 3-inch with 4 adj. dividers per drawer (7.6 cm)
- (2) 6-inch with 4 adj. dividers per drawer (22.9 cm)

#### **Internal Drawer Dimensions:**

• 23" W x 17" D (58.42 cm W x 43.18 cm D)

# **Overall Cart Dimensions:**

- 43" H x 38" W x 21" D (109.2 cm H x 96.5 cm W x 53.3 cm D) w/ pump pole 76" H x 44" W x 21" D (193 cm H x 111.8 cm W x 53.3 cm D)
- Writing surface: 42" H x (106.7 cm)
- Cart net weight: 207 lbs (94.1 kg)
- Approximate shipping weight: 252 lbs (114.5 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D (144.8 cm H x 119.4 cm W x 68.6 cm D)