

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 \* Ph: (719)637-0300 Fax:(719)597-8273 \* www.harloff.com

# Classic Line, Six Drawer, Emergency Cart, Standard Package #6400



Color Shown is Red. Additional Colors Available

## **Specifications**

#### **Construction:**

• Sturdy, uni-body steel construction

#### Features:

- · Painted steel break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Painted steel top rail
- Stabilizing pontoon bumpers
- Five-inch (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- Antimicrobial plastic top, removable
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- · Cart is delivered assembled
- Twelve year warranty on manufactured parts

### **Drawer Configuration:**

- (4) 3"
- (1) 6"
- (1) 12"

#### **Internal Drawer Dimensions:**

23" W x 17" D (58.4 cm W x 43.2 cm D)

#### **Assembled Cart Dimensions:**

- 45.75" H x 32" W x 22" D (116.2 cm H x 81.3 cm W x 55.9 cm D)
- Writing surface: 41.25" (104.8cm)
- Cart weight: 128 lbs (58.6 kg)
- Approx. shipping weight: 180 lbs (81.6 kg)
- Approx. shipping dims: 57" H x 33" W x 26" D (144.8cm H x 83.8cm W x 66cm D)