

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Classic Line, Six Drawer Workstation (Tall) with Key Lock, Standard Package #6450



Color Shown is Beige.
Additional Colors Available

Specifications

Construction:

• Sturdy, uni-body steel construction

Features:

- · Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel one directional, one with brake
- · Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Two aluminum side rails / push handles
- Stabilizing pontoon bumpers
- Twelve year warranty on manufactured parts
- Latex free
- · Durable, hard-baked powder coat finish

Drawer Configuration:

- (4) 3"
- (1) 6"
- (1) 12"

Internal Drawer Dimensions:

• 23" W x 17" D (58.42 cm W x 43.18 cm D)

Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D (107.95 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 129 lbs (58.5 kg)
- Approximate shipping weight: 170 lbs (77.1 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D (134.6 cm H x 83.8 cm W x 66 cm D)