

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

Workstation Cart/ Model: 6451

Six Drawer Workstation with Key Lock, Standard Package



Colors:

Shown in Beige, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

Construction:

Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- · Raised back rail system with three horizontal rails
- Four aluminum side rails / push handles
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel one tracking, one with brake
- · Antimicrobial plastic top, removable
- Stabilizing pontoon bumpers
- · Latex free
- · Durable, hard-baked powder coat finish

Drawer Configuration:

- (4) 3-inch (7.6 cm)
- (1) 6-inch (15.2 cm)
- (1) 12-inch (30.5 cm)

Internal Drawer Dimensions:

• 20" W x 17" D (50.8 cm W x 43.2 cm D)

Overall Cart Dimensions:

- 41.5" H x 32" W x 22" D (105.4 cm H x 81.3 cm W x 55.9 cm D) With Back Rail: 57.75" H (146.7 cm H)
- Writing surface: 40.25" H (102.2 cm)
- Cart net weight: 145 lbs (65.7 kg)
- Approximate shipping weight: 165 lbs (74.8 kg)
- Approximate shipping dims: 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts