

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

6 Drawer Workstation (Tall) with Basic Electronic Pushbutton lock,

Deluxe Package #7456E



Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D (108 cm H x 81.3 cm W x 55.9 cm D)
- With Accessories: 57" H (144.8 cm H)
- Cart net weight: 160 lbs (72.6 kg)
- Approximate shipping weight: 201 lbs (83 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D (134.6 cm H x 83.8 cm W x 66 cm D)

Specifications

Construction:

• Sturdy, uni-body steel construction

Features:

- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel one tracking, one with brake
- · Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Stabilizing pontoon bumpers
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

Specialty Package Accessories Included:

- Raised back rail system with two horizontal rails (68163)
- Four aluminum side rails / push handles (68161)
- 6-Compartment tilt bin organizer (683805)
- Tape and label dispenser (680415)
- Catheter box, 18", with rail clip (680413)
- Large utility basket (680407)
- Plastic waste container, 3 gallon, with bracket and rail clip (684801, 683801)
- Large sharps container and bracket with rail clip (684802, 683801)
- 3-inch exchange trays (2) (681501)
- 6-inch exchange trays (2) (681502)
- 12-inch Adjustable drawer divider set (680503)

Drawer Configuration:

- (4) 3" (7.6 cm)
- (1) 6" (15.2 cm)
- (1) 12" (30.5 cm)

Internal Drawer Dimensions:

• 23" W x 17" D (58.42 cm W x 43.18 cm D)