

OptimAL Line, Aluminum 460 Punch Card Medication Cart #AL460PC

OptimAL Line medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.



Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted aluminum cabinet and drawers

Features:

- Dual column medication cart with 480 punch card capacity, single width narcotics drawer
- BEST® brand lock on cabinet;
- Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Twelve year manufacturer's warranty
- Protective corner bumpers

Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 3 ea - 9.75" punch card drawers

Right Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D (116.2 cm x 91.4 cm x 59.7 cm)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 137 pounds (62.3 kg)
- Approx. shipping weight: 197 pounds (90.9 kg)
- Approx. shipping dims: 57" H x 39" W x 26" D (144.8cm H x 99.1cm W x 66cm D)