

OptimAL Line 600/ Model: AL600PC

Aluminum Punch Card Medication Cart, Standard Package

OptimAL Line medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.



Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

Construction:

- Painted aluminum cabinet and drawers

Features:

- Dual column medication cart with 600 punch card capacity, double width narcotics drawer
- BEST® brand lock on cabinet
- Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Left Column Drawer Configuration:

- (1) - 4-inch multipurpose storage drawer
- (3) - 9.75-inch punch card drawers

Right Column Drawer Configuration:

- (1) - 4-inch multipurpose storage drawer
- (1) - 9.75-inch locking narcotics drawer
- (2) - 9.75-inch multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 111.1 cm W x 61 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 145 lbs (65.9 kg)
- Approx. shipping weight: 190 lbs (86.36 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 60.9 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts