

## OptimAL Line 600/ Model: AL602PC

### Aluminum Punch Card Medication Cart, Standard Package, Double Wide Narcotics Drawer

OptimAL Line medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.



### Construction:

- Painted aluminum cabinet and drawers

### Features:

- Dual column medication cart with 600 punch card capacity, double width narcotics drawer
- BEST® brand lock on cabinet, 4 keys
- Illinois lock on narcotics drawer, 4 keys
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

### Left Column Drawer Configuration:

- (1) - 4-inch multipurpose storage drawer
- (3) - 9.75-inch punch card drawers

### Right Column Drawer Configuration:

- (1) - 4-inch multipurpose storage drawer
- (1) - 9.75-inch locking narcotics drawer
- (2) - 9.75-inch multipurpose storage drawer

### Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 109.2 cm W x 59.7 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 145 pounds (65.9 kg)
- Approx. shipping weight: 190 pounds (86.36 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 60.9 cm D)

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

### Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>